



Validation of anti-KAP1 antibody: rabbit rabbit α -KAP1 (ab10483, Lot # 193224)

- A. Western blot of 10 ug nuclear extracts prepared from the indicated cell lines.
- B. Immunoprecipitation of KAP1 from 293 nuclear extracts. IP-WB of KAP1 Ips using 1, 2, 4 ug of IgG control or anti-KAP1 antibodies. SN is the supernatant from the KAP1 IP reactions. Equal amounts of nuclear extracts were added to each IP reaction.